


ARU[®]


Chemical Products



ruckgaber
orthopädie
service



Contents

	Page
ARU® adhesives	4 - 5
ARU® solvents, last foams	6
ARU® release agents Lupol last foam system	7
ARU® accessories for last foam	8 - 10
ARU® casting resins, fillers, curing agents	10
ARU® fillers, curing agents, liquid toe caps, leather softeners, finishes	11
Additional chemical products from	12 - 13
	
Equipment for working with chemical products	14 - 15

ARU® adhesives

Article	Description	Application	Article no.	Size
ARU-Orthoplast	plastic glue, toluene-free, reactivable, fast reprocessing, good adhesion	pressure: 1-5 bar, exposure time: approx. 5-30 min storage stability: min. 12 months at 20 °C	10633006	0,64 kg
			10633040	4,0 kg
ARU-Orthocolle	polychloroprene glue, transparent, toluene-free, particularly suitable for allergics (free of tert.but. Phenol) and women	also for light-coloured substances, pot life: up to 8 h depending on hardness requirement, pressure: 1-5 bar, exposure time: 10-60 min storage stability: min. 12 months at 20 °C	10634010	0,64 kg
			10634040	4,0 kg
			10634080	8,0 kg
			10634230	23,0 kg
ARU-Ortho TR	polychloroprene glue, transparent, toluene-free, reactivable, fast reprocessing, very good TR adhesion	also for light-coloured substances, glues everything except PVC pot life: up to 4 h depending on hardness requirement, pressure: 1-5 bar, exposure time: approx. 20 min storage stability: min. 12 months at 20 °C	10635040	4,0 kg
			10635080	8,0 kg
			10635230	23,0 kg
ARU-house brand special 170	polyurethane glue, immediate reprocessing, high adhesion, requires activation	specially for absorbent materials, pressure: 1-5 bar, depending on material hardness, exposure time: approx. 20 min storage stability: min. 6 months at 20 °C	10637010	1,0 kg
			10637050	5,0 kg
			10637150	15,0 kg
ARU-rubber solution	natural-rubber base	for all sewing and casting work on shaft and shoes. pressure: 1-5 bar or press on manually exposure time: 10-15 min storage stability: min. 12 months at 20 °C	10638010	0,55 kg
			10638040	4,0 kg
ARU-Craft paste L	glue based on modified dextrin, paste form, high solid content, long open time, does not quickly shrink while drying, very yielding, does not ferment or spoil, high adhesion strength	storage stability: min. 12 months at 20 °C	10669106	0,60 kg
			10669160	6,0 kg



ARU® adhesives

Article	Description	Application	Article no.	Size
ARU superglue 1 , red	Special glue based on cyanacrylate, low viscosity	small repairs of all types, pressure: manual 10 s	10628109	20 gr.
			10628110	50 gr.
storage stability: min. 6 months at 20 °C min. 12 months at 3-6 °C				
ARU superglue 2 , green	Special glue based on cyanacrylate, medium viscosity	small repairs of all types, pressure: manual 10 s	10628107	20 gr.
			10628108	50 gr.
storage stability: min. 6 months at 20 °C min. 12 months at 3-6 °C				
ARU blitz glue 3 , blue	special glue based on cyanacrylate, high viscosity	small repairs of all types, pressure: manual 10 s	10628105	20 gr.
			10628106	50 gr.
storage stability: min. 6 months at 20 °C min. 12 months at 3-6 °C				



overview of adhesives

- dilute with diluter as specified by manufacturer
- dilute with solvent as specified by manufacturer
- dilute with water
- do not dilute
- glue without pretreatment
- glue after pretreatment as specified by manufacturer
- gluing not possible

glue	dilute	leather	rubber	cork	TR	felt	PUR	EVA	PVC	Poro	natural-crepe	textile	foam	Card-board, paper
ARU-Orthoplast	■	●	●	●	●	●	●	●	●	●	●	●	●	●
ARU-Orthocolle	■	●	●	●	●	●	●	●	—	●	●	●	—	●
ARU Ortho TR	■	●	●	●	●	●	—	●	—	●	●	●	—	●
ARU-house brand special 170	■	●	●	●	●	●	●	—	●	●	●	●	—	●
ARU-rubber solution	■	●	—	—	—	●	—	—	—	—	—	●	—	●
ARU-Superglue	■	●	●	●	●	●	●	●	●	●	●	●	—	●
ARU-Craft paste	■	●	—	●	—	●	—	—	—	—	—	●	●	●

The following are available for pretreatment:

Renia RehaGol	KÖ-reinforcer TR 280
10649303 85 ml	10628001 100 ml
10649304 250 ml	10628005 900 ml
Renia curing agent C	KÖ-Pur
10649302 50 ml	10617020 100 ml
10649310 1000 ml	10617022 0,330 kg
Renia Primer	KÖ-Halosol S
10649301 85 ml	10617007 85 ml
10649300 250 ml	10617008 250 ml
	KÖ-Halosol W 5 FL
	10617010 1000 ml



ARU® Solvents

Article	Description	Application	Article no.	Size
ARU-Acetone	cleaning agent	cleaning, pretreatment, greasing, evaporates immediately	10606010 10606050	1,0 l 6,0 l
storage stability: unlimited at 20 °c				
ARU-sole remover	solvent mixture, very penetrating, contains toluene	remove shoe soles, clean heels, clean leather	10639010 10639050	1,0 l 5,0 l
storage stability: unlimited at 20 °c				

ARU® last foams



Article	Description	Application	Article no.	Size
ARU-OP 500 comp. A	PUR last foam, 2-component, impact-resistant, easy to clamp, tack and screw, manual mixing, very easily ground, easy to shape without problems density: 500g/cm ³	manufacture of shoe lasts, orthopaedic lifts, shoe and positive forms, mixing ratio: a:b = 1:1, processing temperature: 18-20 °c	10900015 10900017 10900019 10900021	2,5 kg 5,0 kg 10,0 kg 25,0 kg
storage stability: 6-8 months at 20 °c				
ARU-OP 700 comp. A	see ARU-OP 500, density: 700g/cm ³	see ARU-OP 500	10900025 10900027 10900029 10900031	2,5 kg 5,0 kg 10,0 kg 25,0 kg
ARU-OP comp. B	curing agent for comp. ARU-Op 500 and ARU-OP 700	see aru-op 500	10900026 10900028 10900030 10900032	2,5 kg 5,0 kg 10,0 kg 25,0 kg
storage stability: 6-8 months at 20 °c				



Brafapol L	last cast resin, 2-component, dense structure, easy to clamp, tack and screw, manual mixing, very easily ground, shock, impact and breakage-resistant	manufacture of shoe lasts, orthopaedic lifts, shoe and positive forms, mixing ratio: brafapol: curing powder = 100:8 parts by weight, processing temperature: approx. 20 °c	10600025 10600100	2,5 kg 10,0 kg
storage stability: 6-8 months at 20 °c				
curing powder	for Brafapol L		10600101 10600102	0,2 kg 0,8 kg

ARU® release agents

Article	Description	Application	Article no.	Size
ARU-release agent	wax-based, 1-component, white, cures by evaporation of solvent, economical, simple and easy to process, forms a thin wax layer as release coating	for all hard and soft foam systems, for completely dry plaster, synthetic binders and socks	10900041	1,0 l
			10900045	5,0 l



LUPOL Lupol last foam system

Article	Description	Application	Article no.	Size
Lupol hard foam comp. A	hard-foam system with no foaming factor, 2-component, easy to use, low specific gravity, 3 components with addition of Lupofill	fabrication of positive shapes for lasts and footbeds, e.g. forefoot foambacking, exposure time: approx. 4.5 min, demoulding: after approx. 25 min, fully hardened after 12 h	10901001	1,0 kg
			10901003	2,5 kg
			10901005	5,0 kg
			10901010	10,0 kg
			10901050	50,0 kg



Lupol hard foam comp. B	curing agent for comp. A	exposure time: approx. 4.5 min, demoulding: after approx. 25 min, fully hardened after 12 h, mixing ratio A:B = 1:1	10901002	1,0 kg
			10901004	2,5 kg
			10901006	5,0 kg
			10901011	10,0 kg
			10901051	50,0 kg

Lupofull filler	mineral-based filler, easy to use, low specific gravity, density: 0,05 g/cm³	filler for Lupol hard foam, mixing ratio by weight: 1 part Lupol A, 1 part Lupol B, 1 part Lupofill	10901105	5,0 kg
			10901120	20,0 kg



Lupofull Light filler	filler, low alkali, easy to use, low specific gravity, density: 0,01 g/cm³	filler for Lupol hard foam, mixing ratio by weight: 1 part Lupol A, 1 part Lupol B, 0.2 parts Lupofill Light	10901202	2,0 kg
			10901205	5,0 kg



Luposlip	methanol-based separating agent, very fluid, perfect separating properties, requires Lupotex as base	seals and regulates plaster surface of negative cast	10901301	1,0 kg
			10901302	2,5 kg
			10901305	5,0 kg
			10901310	10,0 kg




Lupotex	latex-based separating agent, very fluid, perfect separating properties	seals and regulates plaster surface of negative cast	10901401	1,0 kg
			10901402	2,5 kg
			10901405	5,0 kg
			10901410	10,0 kg



ARU® accessories for last foam

Article	Description	Application	Article no.	Size
ARU-wood spatulas	short and long models	for stirring cast resins and fillers	10900125	250 mm
			10900130	300 mm




Biplatrix-plaster bandage 	quickset plaster bandage, fast setting, fast reprocessing, sealed in moisture-proof bags, length: 3 m	for negative casts for last and inlay fabrication, immerse in water at 20 °C and press lightly	40000210	10 cm
			40000212	12 cm
			40000215	15 cm




ARU-copy bandage 	synthetic bandage, fast, clean impressions, no folds, castable, length: 3.6 m	fabrication of negatives for last and inlay fabrication, process like plaster bandage, wear safety gloves, avoid direct contact with skin	40000510	10 cm
---	---	---	----------	-------



STS-Quick-Socks 	Quick-Socks OT size 1	cylinder shape	63400061	size S
	Quick-Socks OT size 2		63400062	size M
	Quick-Socks OT size 3		63400063	size L
	Quick-Socks OT size 6		63400064	size XL
	Socklets size 1	short leg	63400071	size S
	Socklets size 2		63400072	size M
	Socklets size 3		63400073	size L
	Socklets size 6	63400074	size XL	
	Quick-Socks foot shape size 1	long leg	63400081	size S
	Quick-Socks foot shape size 2		63400082	size M
Quick-Socks foot shape size 3	63400083		size L	
Quick-Socks foot shape size 6	63400084		size XL	



Foaming aids 	for lasts. 2-part, plastic	for fabrication of forefoot tips	49999971	
--	----------------------------	----------------------------------	----------	--

ARU-mixing beaker 	container for stirring cast resins and foams	cardboard, paraffin-impregnated	55003000	0,5 l
		plastic	55005000	1,0 l
		plastic with scale	55001000	0,4 l

ARU® accessories for last foam

Article	Description	Application	Article no.	Size
---------	-------------	-------------	-------------	------

ARU-latex gloves



disposable gloves,
unsterile, talcumed,
100 units per pack,
sizes:
small, medium, large

57272101	size S
57272102	size M
57272103	size L



ARU-PVC-safety gloves



poly disposable gloves,
unsterile,
100 units per pack,
large

56010100	Size L
----------	--------



Talcum	powder 6/0	10689200	0,8 kg
--------	------------	----------	--------



Talcum dispenser	46288000
------------------	----------



Pevalin Special hand cleaning cream	56005001	200 ml
	56005004	500 ml
	56005005	3,0 l
	56005050	2,2 l



Pevamat dispenser UKS-2 dispenser for 2.2 l cartridge, wall mounting	56005052
---	----------



ARU® accessories for last foam



Article	Description	Application	Article no.	Size
ARU-stretch foil	transparent, very elastic, 30 cm wide, 300 m long	for footbed foam	10980170	
ARU-foil cutter	for stretch foil, wall and table mounting		10980171	
ARU-cleaning paper	roll, 3-ply, blue, 39 cm wide, 1000 sheets 38.5 cm per sheet		10980181	
ARU-roll holder	for cleaning paper roll, wall and table mounting		10980180	

ARU® casting resins, fillers, curing agents



Article	Description	Application	Article no.	Size
ARUfix	polyester casting resin, 2-component, very easy to use	joints, sole reinforcements, reinforcements, processing time: approx. 30 min after adding curing agent	10900070	0,8 kg
curing paste	for ARUfix, white		10900077	25 gr.
Lupokitt	polyester filler, flexible, high adhesion strength, can be nailed	filling corrections on positive forms of hart foam, cosmetic applications, curing additive 1-2%	10901501 10901502	1,0 kg 2,5 kg
curing paste	for Lupokitt, red		51024005	50 gr.
ARU-fibre glass filler	filler with fibre glass reinforcement, polyester resin base, fast curing, high resistance to water, petrol, oils, acids and alkalis, very high strength, easy to grind, green-yellowish-transparent	joint reinforcement, installation of joint springs, sole reinforcements, supports and reinforcements, curing additive approx. 2%, processing time: approx. 6 min after adding curing agent, grind after 20-30 min	51072415	1,5 kg
curing paste	for ARU fibre glass filler, red		51024005	50 gr.

ARU® fillers, curing agents, liquid toe caps, leather softeners, finishes

Article	Description	Application	Article no.	Size
ARU-light filler SL	polyester-based filler mass, extremely light, easily vacuum-formed and ground, low dust during grinding, high productivity, pore-free, very good adhesion, can be nailed and clamped, beige	filling corrections on positive forms and hard foam, cosmetic applications, curing agent: approx. 2%, processing time: approx. 6 min after adding curing agent, grindable after 15-25 min	51072310	1,0 kg
			51072313	1,6 kg
			51072317	7,5 kg
curing paste	for ARU light filler SL, red		51024005	50 gr.
curing paste	for ARU light filler cartouche, red		51024004	40 gr.
ARU-dispenser	for sl filler mass, mobile or fixed installation, for clean dosage, easy to use, empties completely, curing agent added automatically at 2%, complete with curing paste and filler		10980160	
ARU-liquid shoe caps	easily applied	for thinly coated shoe caps	10900051	1,0 kg
			10900052	5,0 kg
			10900053	15,0 kg
ARU-leather softener	makes all walking, stretching, moulding and gusseting processes easy	for all leather	10900101	1,0 l
			10900105	5,0 l
Top Finish	edge and sole finishing dye, good coverage, silky gloss, water-soluble, mixable, dries fast, fast colours when dry	colours: black dark brown colourless	10674001	0,5 l
			10674005	0,5 l
			10674022	0,5 l
Top Finish pad	sponge applicator with valve, makes application easier	colours: black 081 medium brown 035 dark brown 046 dark blue 078 mud 044 nickel 011 grey beige 017 bison 022 colourless 100 light beige 05 white 09 stone grey 056 light grey 060 sand 07 smoke 054	10674101	0,75 ml
			10674104	0,75 ml
			10674105	0,75 ml
			10674107	0,75 ml
			10674115	0,75 ml
			10674116	0,75 ml
			10674117	0,75 ml
			10674120	0,75 ml
			10674122	0,75 ml
			10674129	0,75 ml
			10674130	0,75 ml
			10674131	0,75 ml
			10674137	0,75 ml
			10674138	0,75 ml
10674139	0,75 ml			
replacement applicator	for round bot		10674199	

cartouche 51072310





COLLE DE COLOGNE NEU

10640006	0,64KG
10640040	4 KG
10640230	23 KG

RENIA TOP FIT

10641001	60 GR.
10641006	0,64 KG
10641040	4 KG
10641080	8 KG

RENIA FORTE

10642080	8 KG
----------	------

RENIA VULKOFEST 96

10643006	0,64 KG
10643040	4 KG
10643080	8 KG
10643230	23 KG

RENIA MULTICOLLE

10644006	0,64 KG
10644040	4 KG

RENIA NOVO COLLE

10645080	8 KG
----------	------

RENIA RENOL -ASTI-

10646006	0,64 KG
10646040	4 KG

RENIA SYNTIK TOTAL

10647001	90 GR
10647010	1 KG
10647040	4 KG

RENIA AQUILIM

10648005	0,5 KG
10648030	3 KG

RENIA QUICK SOHL TUBE

10649001	black/brown
10649002	leather colour
10649003	old white

RUBBER SOLUTION -EXTRA-

10649006	0,6 KG
----------	--------

RENIA PRIMER

10649301	85 ML
10649300	250 ML

RENIA CURING AGENT -C-

10649302	50 ML
10649310	1000 ML

RENIA REHAGOL

10649303	85 ML
10649304	250ML

DILUTER f. FAST WORKING

10649406	
10649450	

DESOHL SOLE REMOVER

10649510	1 LTR
10649550	5 LTR

RUBBER SOLUTION

10649650	4 KG
10649653	23 KG

EDELKITT-LEATHER GLUE

10649750	4 KG
----------	------

RENIA ORTEC SPECIAL

10650006	0,64 KG
10650040	4 KG
10650230	23 KG

RENIA ORTEC SPECIAL

10650006	0,64 KG
10650040	4 KG
10650230	23 KG

HOOCO FINISHING WAX

10686001	black
10686003	light brown
10686004	medium brown
10686005	dark brown
10686019	red
10686020	blue
10686022	colourless
10686026	negro
10686027	red brown
10686500	optima

RENIA LEATHER FINISH

10688001	50 ML
10688003	100 ML
10688010	1 LTR

BRILLANT FINISHING WAX

10689001	black
10689004	medium brown
10689005	dark brown
10689022	colourless
10989026	negro

KÖ-SPECIAL-ADHESIVE

10611006	0,64KG
10611040	4 KG

KÖ-MEGABOND

10612040	4KG
10612100	10 KG
10612230	23 KG

KÖ-HALOSOLS

10617007	85 ML
10617008	250 ML

KÖ-HALOSOL W 5 FL

10617010	1 LTR
----------	-------

KÖ-PUR UNDERCOAT

10617020	100 ML
10617022	330 GR

KÖ-ANTISTATIC FL

10617025	250 ML
----------	--------

KÖRABOND 11

10617030	100 ML
----------	--------

KÖ LIGHT FILLER

10617050	0,5 KG
----------	--------

CURING AGENT FOR LIGHT FILLER

10617051	
----------	--

KÖRACRAFT

10618040	4 KG
----------	------

KÖVULFIX REKORD

10620006	0,64 KG
10620040	4 KG

KLEBSTOFF KÖ F 6

10621040	4 KG
----------	------

KÖRAPLAST SF

10626006	0,64 KG
10626020	2 KG
10626040	4 KG
10626100	10 KG

KÖRASOLF SOLVENT

10627001	PU 1 KG
10627002	CR 1 KG
10627050	PU 5 KG
10627055	CA 5 KG

KÖ REINFORCER TR 280

10628001	100 GR
10628005	800 GR

KÖ DILUTER

10628006	1 LTR
10628050	5 LTR

KÖ-ATOMGLU

10628095 C 2	yellow 20 GR
10628130 C 3	blue 20 GR
10628104 C 3	blue 50 GR

HOOCO NEOTEX EXTRA

10630040	4 KG
10630230	23 KG

HOOCOFIX THERMO CONTACT

10631040	4 KG
----------	------

BKF LAMINATING 80:20

51014010	900 GR
51014046	4,6 KG

BKF POLYFEST 98

51015010	1 KG
51015030	5 KG

BKF KNET SILICON

51016030	Komp.A 300 GR
51016031	Komp.B 300 GR

ASTI OP

51017050	Komp.A 4,6 KG
51017051	Komp.B 4,6 KG

CURING PASTE RED

51024005	50 GR
51024020	200 GR
51024050	500 GR

CURING POWDER

51025030	150 GR
51025040	500 GR

ARALDIT CASTING RESIN OP

51042010	1 KG
----------	------

CURING AGENT HY FOR

ARALDITE

51042040	HY 956 200 GR
51042080	HY 943 200 GR

LAST CASTING COMPOUND

BKF A

51045000	5 KG
----------	------

CLEANING AGENT W

51045001	1 LTR
----------	-------

AEROSIL THIXOTROPIC AGENT

51072042	500 GR
----------	--------

FILLER COMPOUND C

51072231	1,0 KG
51072232	2,5 KG
51072235	5,0 KG

HARD FOAM V

52002011	A 0,910 KG
52002012	B 1,090 KG

ASTI TECH

52017011	soft A 0,865 KG
52017051	hard A 0,865 KG
52017091	cur.agentB0,865 KG

ANTICAKING AGENT 65/98

52023000	1 LTR
52026050	5 LTR

ANTICAKING AGENT 65/98

52031700	150 GR
52031707	750 GR

WETTING AGENT FOR

ANTICAKING AGENT 65

52031810	
----------	--

ANTICAKING AGENT 69

52034010	1 LTR
52034020	2 LTR

ANTICAKING AGENT MO 70

52034071	1 LTR
52034095	5 LTR



FATO INTERIOR GLUE 142

61100142 1 KG
61100145 5 KG

FATO HIGH-GLOSS INK

61201001 black

FATO SILKY LUSTER

61201001 black
61203005 medium brown

FATO PVC DYE

61204001 black

FATO FINISHING DYE

61209001 black
61209004 medium brown
61209005 dark brown
61209008 havanna
61209013 mahogany

FATO TRANSPARENT FINISHER

61209010 natural
61209011 teak
61209012 rosewood

FATO SPRAY POLISH

61241501 black 5 LTR
61241510 natur 5LTR

FATO REPAIR PENS

61251001 / unit
61251012 / dozen
61251020 / wood box

FATO GRINDING WAX

61270001 black
61270005 dark brown
61270008 havanna
61270010 natural

FATO POLISHING WAX

61271010 Colourless

FLEX LEATHER STRETCHER

61281050 6 LTR

STALIT KG

63100010

CURING AGENT FOR STALIT/MERCALITE

63101000 red 50 GR

TUCA 1402 ANTICAKING AGENT

63140010 1 LTR
63140050 5 LTR

SR LEATHER SOFTENER

63202010 0,75 LTR
63202050 5 LTR
63202300 30 LTR

SILONDA AN 4

63812000 500 ML

SPITACID AS 24

63821000

INCIDIN M-SPRAY EXTRA YMA 24

63822000

MANIPUR AD 3

63825020 6 LTR



EURO ATOMGLU

10628101 50 ML
10628102 blue 50 ML

FULLFIX

10629890 500 GR

HALOPRIM

10657002 250 ML
10657005 500 ML

PRENOFIX 80/10

10661040 4 KG
10661080 8 KG
10661230 23 KG

PRENOFIX UNIPLAST 33

10662040 4 KG

PRENOFIX 100

10664040 4 KG
10664080 8 KG
10664220 22 KG

VA VITE SOLVENT

10665006 1,0 LTR
10665035 3,5 LTR

EUROPREN 2618

10667205 5 KG
10667210 10 KG



LEISTON NEW

57825205 5 KG

CURING AGENT FOR LEISTON

57825215 150 GR

FILLER FOR LEISTON

57825275 750 GR



COLD POLISHING INK

10670301 black 1 LTR

FINISHER 2000

10673101 black 100 ML
10673105 d'brown 100 ML
10673125 optima 100 ML

SEWING MACHINE OIL

10678005 0,5 LTR

MORELLO EXTRA BLACKENER

10681200 50 ML
10681301 1 LTR

BEE SWAX

10684000 RIEGEL

MORELLO LEATHER OIL

10685005 0,5 LTR

MORELLO LEATHER DYE

10687501 black 40 ML
10687504 m'brown 40 ML
10687505 d'brown 40 ML
10687506 beige 40 ML
10687520 blue 40 ML
10687524 d'blue 40 ML
10687525 yellow 40 ML
10687526 grey 40 ML
10687529 linen 40 ML
10687530 white 40 ML
10687555 red 40 ML



Working with chemical products such as glues, casting resins and foams involves risk to health. This applies both for processing and correct storage.

Optimum ventilation in the workplace and storage of hazardous substances in accordance with the regulations are the essential prerequisites for modern, rational and healthy work practices in all workshops.

We offer the right equipment for all areas. We support you with planning and will set up your workplace as specified by you. Turnkey solutions are also available.

Please contact us to make an appointment.

Safety cabinet

manufactured in conformity with DIN 12925 (T1) and TRbF 20 (System L) with self-actuating fire doors, integrated floors with high load capacity, available in various sizes.

Resin cabinet

Compact design with integrated fan exhaust motor, working base for resin casting, storage cabinet, roof doors.



binets

sin casting workstation

esign, with
an, filter system,
nted on swivel arm,
e specially coated
ing, material
nets with dustp

Gluing workstations

stable steel plate construction,
with integrated fan, filter system,
multiplex work base,
large exhaust funnel,
workstation height-adjustable,
optionally with shelf,
available in various sizes

Exhaust systems

we supply the complete piping system and accessories





Ruckgaber Orthopädie Service GmbH
Roseckstraße 35
Gewerbegebiet Seebronn
D-72108 Rottenburg

Tel.: +49(0)74 57 / 94 97-0
Fax: +49(0)74 57 / 94 97-35

E-mail: info@ruckgaber.de
www.ruckgaber.de